

EL DORADO COUNTY EMS AGENCY

FIELD PROCEDURES

808

Effective: July 1, 2011
Last Revised: August 2024
Scope: ALS

(on file)

EMS Agency Medical Director

ETCO2 MONITORING

PURPOSE

The measurement of End-Tidal CO₂ (EtCO₂) is the optimal method of continuously monitoring systemic perfusion in both adult and pediatric patients. It measures expired carbon dioxide using infrared spectroscopy.

INDICATIONS:

Initial and continuous EtCO₂ waveform monitoring must be employed on all patients with advanced airway interventions including:

- Endotracheal Tube
- Supraglottic Airway (SGA)
- Needle Cricothyrotomy
- Respiratory Distress/Arrest

CONTRAINDICATIONS:

- None

PROCEDURE:

1. Select appropriate sized sensor based on the size/type of airway device selected.
2. Calibrate with room air prior to placing on device. Monitor will show a reading of "0"
3. Place sensor between the end of the advanced airway device and the bag valve in accordance with manufacturer specifications.
4. Ventilate patient to maintain a normal ETCO₂ waveform and numeric value (Target = 35-45 mmHg).

- a. Certain conditions that cause metabolic acidosis or poor perfusion, such as hemorrhage, sepsis, DKA, etc., will show lower values even when ventilation is adequate.
5. Generally a normal ETCO₂ waveform (square shaped) should be noted for confirmation of ETT placement. This will typically correspond with a value of 35-45 mmHG or 10-20 mmHg during effective CPR.
6. Changes in waveform may reflect the complications below:
- a) Tube displacement or esophageal intubation – diminished or no waveform
 - b) Obstructed airway or ventilation device failure – diminished or no waveform
 - c) Hyperventilation – low numeric value
 - d) Hypoventilation – high numeric value
7. Include ETCO₂ waveforms from transport and transfer of care in the ePCR

SPECIAL CONSIDERATIONS:

- A drop in ETCO₂ during CPR may indicate suboptimal compressions.
- A sudden rise in ETCO₂ during CPR may be an indicator of "Return of Spontaneous Circulation" (ROSC) and the patient should be re-evaluated at that time.
- Low ETCO₂ levels may be noted in states of profound shock. Ventilation changes should not be done to try and correct ETCO₂ levels but rather correct the underlying cause of the hypo-perfusion state.